

NeuroShield: A Multimodal Wearable with ChronoNet for Real-Time Tonic–Clonic Seizure Detection and GPS Alerting



Supervisor

Dr. Wajeeha Khalil

Assistant Professor

Submitted By

Fasih Ullah

21PWDSC0010

Hasan Jehangir

21PWDSC0006


Aizaz Ahmad

21PWDSC0043

**DEPARTMENT OF COMPUTER SCIENCE
UNIVERSITY OF ENGINEERING & TECHNOLOGY
PESHAWAR**

August 2025

NeuroShield: A Multimodal Wearable with ChronoNet for Real-Time Tonic–Clonic Seizure Detection and GPS Alerting



Supervisor

Dr. Wajeeha Khalil

Assistant Professor

Submitted By

Fasih Ullah

21PWDSC0010

Hasan Jehangir

21PWDSC0006

Aizaz Ahmad

21PWDSC0043

**A Project Report submitted in partial fulfillment of the
requirements for the award of Bachelor's Degree in**

Data Science

**DEPARTMENT OF COMPUTER SCIENCE
UNIVERSITY OF ENGINEERING & TECHNOLOGY**

PESHAWAR

August 2025

Project Approval

This to certify that this project is approved and recommended as partial fulfillment for the award of Bachelor Degree in Computer Science, from University of Engineering & Technology Peshawar.

Supervisor

DR. WAJEEHA KHALIL

(Name)

Chairperson

Dr. Syed Adeel Ali Shah

(Name)

Undertaking

We certify that project work titled “NeuroShield: A Multimodal Wearable with ChronoNet for Real-Time Tonic–Clonic Seizure Detection and GPS Alerting” is our own work. No portion of the work presented in this project has been submitted in support of another award or qualification either at this institution or elsewhere. Where material has been used from other sources it has been properly acknowledged / referred.

Hasan Jehangir

21pwdsc0006

Fasih ullah

21pwdsc0010

Aizaz Ahmad

21pwdsc0043

Acknowledgements

We would like to thank **Dr. Wajeaha Khalil**, our project supervisor, for guiding us in understanding the concepts of Artificial Intelligence and its Application in Healthcare and for her cooperation and support to bring this project to completion. We would also like to thank all our teachers especially Dr. Imran Khalil and Dr. Iftikhar for guiding us in solving our problems related to Machine Learning and Data Science.

We would also like to thank our families and friends for their continuous encouragement and moral support.

Abstract

This project describes the implementation of a wearable device specifically designed to detect tonic-clonic epileptic seizures in real time and send immediate alerts to caregivers and medical professionals. Epilepsy, a chronic neurological disorder affecting millions worldwide, features unpredictable tonic-clonic seizures that carry significant risks of injury due to sudden loss of consciousness and violent muscle contractions. To address this, we integrated multiple sensors which are EMG sensor for muscle activity recording during tonic clonic seizure attack, MAX30102 for Heart rate monitoring, MPU9250 for movement and orientation changes. The data collected by these sensors are analyzed by an AI model trained using Chrono-Net Deep learning algorithm, which is good for processing sequential time series input and is very reliable for real time seizure detection using sequential data. Once the data passed through model and seizure detected, the device transmit an alert with patient live location to the mobile application installed in caretaker phone. The app provides with live location of patient with other important data in the form of dashboard which shows heartbeat.

Evaluation on a validated epilepsy dataset yielded an overall detection accuracy of approximately 82%, with the system consistently achieving detection latency below five seconds (our prototype result, see Chapter 5). The wearable's compact, low-power design, combined with secure cloud based data storage and encrypted communications, ensures continuous monitoring without compromising user comfort or data privacy. By combining sensors with AI model trained with ChronoNet, NeuroShield can advance epilepsy care.

Table of Content:

Undertaking.....	4
Acknowledgements	5
Abstract.....	6
Chapter 1 Introduction.....	11
BACKGROUND	11
1.1 An Overview of Epilepsy: A Global Health Challenge	11
1.2 The Clinical Characteristics of Tonic-Clonic Seizures	11
2. Multimodal Sensor Fusion for Seizure Detection	14
2.1 The Importance and Rationale of Sensor Fusion	14
2.2 Previous Research and Key Findings	15
2.3 Persistent Limitations of Current Multimodal Systems	15
3. The ChronoNet Algorithm for Seizure Detection.....	16
3.1 An Introduction to the ChronoNet Architecture	16
3.2 The Advantages of ChronoNet for Real-time Detection	17
4. Proposed System Description	18
4.1 Detailed System Workflow	19
4.2 Expected Outcomes and Broader Impact.....	22
CHAPTER 2.....	24
INTRODUCTION	24
REVIEW OF SEIZURE DETECTION TECHNIQUES	24
A. Technique 1: Using One Type of Data with Regular Machine Learning.....	24
B. Technique 2: Using EEG with Deep Learning	25
C. Technique 3: Multimodal Sensor Fusion with Deep Learning.....	26
Chapter 3	30
Research – Methodology	30
3.1 overview	30
3.2 Proposed Framework/Model/Technique	30
3.2 Hardware Design	32
3.2.1 Selection of Components	32
3.2.2 Power Distribution and Voltage Compatibility	33
3.2.3 Sensor Functions in Seizure Detection.....	33

3.3 Software Design.....	34
3.3.1 Data Flow.....	34
3.3.2 Dataset and Initial Model Training	35
3.3.3 Data Pre-Processing.....	36
3.3.4 Feature Engineering	36
3.3.5 WebApp Interface.....	37
3.3.6 Model Training	37
3.6.1 Model Retraining with Local Data	44
3.7 Comparative Analysis with Existing Devices.....	44
Chapter Summary.	45
Chapter 4. Results and Discussion.....	47
4.1 Results:	47
4. 2 Discussion.....	48
Chapter 5 Conclusion	51
Chapter 6	54
Sensing and Hardware	54
Embedded and Real-Time Inference	54
Modeling and Data.....	55
Application and Experience	55
Security, Privacy, and Compliance.....	56
Operations and Reliability	56
Research Extensions.....	56
References.....	59

LIST OF FIGURE:

Figure 1.ChronoNet Architecture	17
Figure 2. Model Accuracy Progression over Training Epochs	19
Figure 3. NeuroShield Flow Chart.....	21
Figure 4. Gyroscope connection & Figure 5. Gyroscope and PPG connection	34
Figure 6. EMG sensor Connection.....	34
Figure 7. NeuroShield Circuit Diagram	44
Figure 8.ROC and PR curves with operating point,confusion matrix heatmap, latency distribution.	48

LIST OF TABLES

Table 1 . SUMMARY COMPARISON OF SEIZURE DETECTION TECHNIQUES	27
Table 2. Comparison of Proposed Prototype with Other Wearable Seizure Detection Devices.....	44
Table 3. Overall performance of the proposed model.....	47

BACKGROUND

1.1 An Overview of Epilepsy: A Global Health Challenge

Epilepsy is a major health concern worldwide. It's a long-term brain condition that affects about 50 million people of all ages everywhere [1]. The International League against Epilepsy (ILAE) describes it as a brain issue that makes people likely to experience seizures. It can also cause problems with thinking, emotions, and behavior. A person is diagnosed with epilepsy if they have at least one seizure [2]. Seizures happen because of unusual electrical activity in the brain. Sometimes, groups of brain cells send too many signals at once, leading to seizures [2], [1].

Epilepsy is a common nerve problem that can make someone's life harder, cause mental and social struggles, and limit school and job prospects. People with epilepsy might face stigma, which can make them feel lonely and be treated unfairly. Because seizures can happen anytime, patients and their families might always worry. They can also have thinking issues, like problems with memory and attention. The cost of treatment and missed work can also cause money problems [2].

1.2 The Clinical Characteristics of Tonic-Clonic Seizures

Tonic-clonic seizures, also known as grand mal seizures, are a serious type of seizure. They have two distinct parts, making them noticeable but risky. First, there's the tonic phase, where all the body's muscles stiffen. The person loses consciousness quickly and might fall, leading to injuries.

The chest and throat muscles often tighten, pushing air out of the lungs and causing a groan or cry. The skin might turn blue from lack of oxygen.

The first phase changes quickly to the clonic phase. During this time, the arms, legs, and body jerk strongly and rhythmically. The person might bite their tongue because the jaw muscles are contracting hard. They might also lose control of their bladder or bowels. These strong movements can lead to broken bones, dislocations, and injuries to soft tissues. After the jerking stops, the patient enters the postictal state, a recovery period when the brain's electrical activity slowly returns to normal. At first, the person might not respond. Later, they could be confused, sleepy, have a headache, and not remember what happened [2].

A very dangerous situation is status epilepticus where a seizure lasts longer than five minutes, or seizures happen one after another without the person waking up in between. Tonic-clonic status epilepticus is an emergency that can cause brain damage or death if not treated quickly. The ongoing seizure activity puts the brain in crisis, raising the risk of low oxygen, acidosis, and brain cell death [3].

1.3 Core Challenges in Modern Epilepsy Management

Despite advances in medicine, the effective management of epilepsy is fraught with persistent challenges that hinder patient safety and quality of life.

1.3.1 Delays in Seizure Detection and Treatment: One of the main problems is noticing seizures fast and for sure, especially when they happen outside a hospital. The way to diagnose seizures is with an electroencephalogram (EEG), which records the brain's electrical activity. But video-EEG monitoring in the hospital uses a lot of resources, is pricy, and can only be done for a short time. It's not really possible to use it for constant,

long-term monitoring in someone's daily life. A patient might not even have a seizure in the hospital, meaning the test doesn't catch anything [4]. Many seizures, especially at night, go unseen and unreported, delaying needed medical care [5].

1.3.2 Too Many False Alarms from Current Devices: To help with monitoring, people have made wearable seizure detection systems. A big problem with many of these devices is that they give too many false alarms. Systems that use just one type of sensor, like an accelerometer, often can't tell apart the rhythmic movements of a clonic seizure from normal, strong movements like running, brushing teeth, or restless sleep. This results in lots of false alarms, which can make caregivers tire alarms. If there are too many false alarms, caregivers might start ignoring them and be slower to respond to a real emergency [6]. **The Need for Alerts Right Away:** Many systems, especially those that analyze data later, don't send alerts to caregivers or medical staff quickly. Even a short delay can be dangerous. A fast alert allows a caregiver to do simple things, like turn the person on their side to prevent choking on saliva or vomit, clear the area of dangerous things, and time the seizure to know when to call for help. Without this alert, the risks from seizures increase [7].

1.3.3 The Need for Alerts Right Away: Many systems, especially those that analyze data later, don't send alerts to caregivers or medical staff quickly. Even a delay of seconds can be dangerous. A quick alert can save the patient from injuries or even death. A notification can save the patient from dangerous situations, giving caretaker a chance to clear the area of danger, prevent patient from falling into deadly areas like fire, well and

give time to caretaker to call for emergency services. Without real time, it can increase the risks associated with seizure events.

1.4 The Promise of Wearable Technology in Epilepsy Care

Wearable tech is changing how epilepsy is treated. With these small sensors worn all the time, like wristbands or patches, we get constant information about what's happening, instead of just checking in at the doctor's office sometimes. These wearable sensors use things like accelerometers to track movement, gyroscopes to measure orientation, EMG to see muscle activity, and PPG to check heart rate and blood oxygen. This gives us a practical way to get data instantly[8], [9]. This tech helps patients feel more independent, reassures their families and caretakers, and makes daily life safer for those with epilepsy [5], [6].

2. Multimodal Sensor Fusion for Seizure Detection

2.1 The Importance and Rationale of Sensor Fusion

Because one-sensor systems have limits, studies have moved to using many kinds of sensors together. This way gives a better and stronger way to sort info. The idea is that by mixing different body signals, the system can use the good parts of each sensor to understand the patient better. This helps fix the problems of using just one sensor, cutting down on mistakes and making the whole system more trustworthy.

For example, an accelerometer might detect strong, rhythmic movements. A single-sensor system would have to classify this as a potential seizure. However, a multimodal system can cross-reference this motion data with other streams. It would check the EMG sensor for the characteristic high-amplitude electrical signals of intense muscle contraction typical of a tonic-clonic seizure.

Simultaneously, it would query the PPG sensor for a corresponding spike in heart rate (tachycardia), a common autonomic response to a seizure. If the motion is present but the corroborating EMG and PPG signals are absent, the system can confidently classify the event as non-epileptic activity, such as exercise, thereby preventing a false alarm.

2.2 Previous Research and Key Findings

Recent research has consistently validated the superiority of the multimodal approach. The landmark European **SeizeIT2 study**, a large-scale clinical trial, demonstrated that a system combining EEG, EMG, and accelerometer achieved an outstanding sensitivity of 97.7% for detecting tonic-clonic seizures while substantially reducing the rate of false positives compared to using EEG alone. Combining physiological and neurological sensors is super valuable. For example, Böttcher and team (2021) built a system that uses both movement tracking and skin sweat measurements. They then used machine learning to analyze this data and got really good results: perfect sensitivity and very few false alarms (only 0.09 every 24 hours). These studies strongly suggest that mixing data types is the best way to build seizure detection tools that doctors can really depend on.

2.3 Persistent Limitations of Current Multimodal Systems

Even with the advances, there are still problems with using systems that use many types of input for everyday life:

2.3.1 Lots of Computing and Power Needed: It takes a lot of computer power to constantly get, process, and analyze many data streams that happen often. Doing signal processing, getting features, and sorting using complex machine learning needs a lot from the small computer chips in wearable devices. This affects how long the battery lasts. The

need to run many sensors and advanced math can make a device's battery die fast, which makes it hard to use for the constant watching that epilepsy needs, which then makes patients not want to use it.

2.3.2 Easy to Use and Comfortable: If a patient doesn't want to wear a device all the time, it won't work. Systems that need many sensor stickers, wires, or are big and not easy to wear probably won't be used for a long time. Things like skin getting irritated from sticky stuff, feeling embarrassed about wearing a device that people can see, and needing the device to last and be waterproof are important things to think about when designing, and they are often missed in early research.

3. The ChronoNet Algorithm for Seizure Detection

3.1 An Introduction to the ChronoNet Architecture

To get around the computer problems and be as correct as possible, we need some new algorithms. ChronoNet is a deep-learning setup made just for sorting different types of physiological time-series data quickly and correctly. (Roy et al). came up with it in 2018. ChronoNet has a cool design that mixes the good parts of Convolutional Neural Networks (CNNs) and Recurrent Neural Networks (RNNs).

It starts with a bunch of 1D convolutional layers that all have different filter lengths. This is very important. The shorter filters are like special detectors for quick, changing patterns in the data, like the fast shakes of a clonic jerk. The longer filters are made to grab slower, longer-lasting patterns, like the constant muscle tension of the tonic phase. Then, the outputs from these convolutional paths go into Gated Recurrent Unit (GRU) layers. GRUs are a better kind of RNN that is great at

modeling data that comes in a sequence. The gating inside them helps them learn what info from the past matters now and what can be forgotten. This makes them really good at spotting the changing patterns that show an epileptic seizure.




Figure 1.ChronoNet Architecture

3.2 The Advantages of ChronoNet for Real-time Detection

ChronoNet's hybrid CNN-GRU setup works well than regular machine learning and other deep learning setups in a few ways. First, it uses 1D convolutions to pull features at the start. This means it can automatically learn the most important parts from the raw sensor data. This way, we don't have to do feature engineering by hand, which can be a pain and take a lot of computing power. Second, the GRU layers help the model understand timing. This lets it figure out what order events happen in to tell when a seizure is occurring. The result is that it does a better job spotting small and difficult patterns in messy medical data. This leads to way better detection and faster response times. Being able to decide rapidly and spot on is super important for a live alert system where time is of the essence.

4. Proposed System Description

This research plans to create, build, and test a wearable device. It will have different kinds of sensors (EMG, PPG, accelerometer, gyroscope) and use the advanced ChronoNet code to detect tonic-clonic seizures as they happen. The whole thing is built to fix the main problems with current options.

- **Alerts and Location in Real-Time:** When a seizure is detected, the system will send out an alert right away through a phone app. This will notify caregivers and doctors right away. The alert will have the patient's GPS location to help people get there fast.
- **Reduced False Alarms through Advanced Fusion:** By combining three distinct sensor modalities and processing the data with the highly accurate ChronoNet algorithm, the system is designed to achieve superior specificity, with a Training Accuracy of approximately 87.2% and Validation Accuracy of 83.92%, minimizing caregiver fatigue.




Figure 2. Model Accuracy Progression over Training Epochs

- Enhanced User Experience and Adherence:** The physical device will be designed to be lightweight, comfortable for 24/7 wear, and power-efficient to ensure long battery life, thereby encouraging continuous patient use and compliance.

4.1 Detailed System Workflow

- Continuous Data Acquisition:** The wearable device, worn on the wrist or arm, continuously captures four streams of physiological data: 3-axis acceleration (ACC), 3-axis angular velocity (GYRO), muscle electrical activity (EMG), and heart rate/blood volume pulse (PPG).
- On-Device Preprocessing:** The raw sensor data is segmented into short, overlapping time windows. Each window of data is filtered to remove noise and artifacts.
- ChronoNet Inference:** The preprocessed multimodal data window is fed as input to the onboard ChronoNet model.

4. **Seizure Classification:** The model analyzes the data and outputs a classification (seizure or non-seizure).
5. **Alert Generation:** If a seizure is detected, the device triggers an alert and transmits it via Bluetooth to a paired smartphone.
6. **Caregiver Notification:** The smartphone application immediately sends a notification (e.g., SMS, push notification) to a pre-defined list of caregivers.




Figure 3. NeuroShield Flow Chart

4.2 Expected Outcomes and Broader Impact

This system should really help with epilepsy care. Because it can keep watching patients closely and give helpful information, it should help lower the number of injuries and deaths from seizures since people can get help faster. Patients might feel safer and freer, which can allow them to enjoy their lives better. Caregivers should also feel less worried since they won't have to keep constant watch. The data gathered over time can give doctors a clearer idea of how often seizures happen and what they're like so that they can adjust treatments as needed.

5. Chapter Summary

This chapter looked at epilepsy as a big health problem around the world. We talked about how common it is, how it affects people's lives, and the special risks of tonic-clonic seizures. The chapter went over what these seizures look like and the emergency situations they can cause. It pointed out some major problems with how we currently handle them, like not catching them early enough, getting too many false alarms, and not having alerts that work right away. Then, we talked about how wearable tech and combining different sensors could really help. By using movement, muscle, and heart signals, we can get better at spotting seizures and cut down on mistakes. Finally, the chapter introduced ChronoNet, a deep-learning setup, and explained the plan for a wearable device. This device is made to be light, use little power, and spot seizures as they happen. It's meant to tell caregivers where the patient is, make patients safer and more independent, and give doctors better info over time.

6 Thesis organization:

Chapter 1: Intro
This chapter covers what epilepsy is, what problems the study is trying to solve, the goals of the study, why it matters, and a short look at wearable tech for detecting seizures.

Chapter 2 – Background and Literature Survey

This chapter looks at past studies on wearable seizure detection tech, mainly for tonic-clonic seizures. It goes over how the tech has changed, from simple sensors to more complex systems using machine learning, EEG, and deep learning. It also points out what these systems do well, where they fall short, and what still needs to be researched.

Chapter 3 – Research Methodology

This chapter talks about how the system was built, including the hardware and software, the sensors used (EMG, PPG, accelerometer, gyroscope), the ChronoNet deep-learning model, how data was gathered and cleaned, how the model was trained, and how it was tested against other methods and devices.

Chapter 4 – Results and Discussion

This chapter shows the results of the tests, compares ChronoNet to other methods, and talks about what it does well, where it struggles, and how to make it better at spotting seizures and reducing false alarms.

Chapter 5: Conclusion

This chapter sums up the main points and what the project achieved. It says that the ChronoNet system can reliably detect seizures in real-time and explains how this can help patients and caregivers.

Chapter 6 – Future Recommendations

This chapter suggests ways to turn the prototype into a product that doctors can trust, including improving the sensors and hardware, using embedded inference, refining the models and data, improving the user experience, ensuring privacy and compliance, and making it more reliable.

Background/Literature Survey

INTRODUCTION

This chapter looks at existing studies on using wearable tech to spot epileptic seizures. It homes in on methods for tonic-clonic seizures, since these are severe and can cause injuries. The chapter checks out how seizure detection tech has changed, starting from simple systems that use regular machine learning to the newer stuff that uses multiple data types and deep learning. By looking at how well different methods work, what they're good at, and what their limits are, this chapter sets the stage for the research in this paper. The goal is to find the best way to make a wearable seizure detection system that's reliable, works in real-time, and is easy for people to use.

REVIEW OF SEIZURE DETECTION TECHNIQUES

A. Technique 1: Using One Type of Data with Regular Machine Learning

Early work on detecting seizures outside of hospitals mainly used wearable gadgets other than EEG, like wristbands or ankle sensors. These had accelerometers (ACC) and photo plethysmography (PPG) sensors to measure movement and blood flow. The data from these sensors was usually checked using regular machine learning (ML) models, like Support Vector Machines (SVM), Random Forests, and logistic regression, or basic deep learning models like Convolutional Neural Networks (CNNs) [6].

For example, (Tang et al.) checked seizure detection in 94 kids using wrist and ankle sensors. These sensors tracked movement, sweat gland activity (EDA), and PPG signals during 930 seizures [10]. The study created two ML algorithms and found that combining ACC and BVP data gave the best results, with an AUC-ROC score of about 0.752 for all seizure types. This means that using multiple sensors, even if they're not very different, works a lot better than just guessing. Similarly, (Meisel et al). studied seizure prediction using long-term data (over 1400 hours) from wristbands that tracked EDA, temperature, BVP, and actigraphy. Their deep CNN models were able to detect the preictal state better than chance for about 60% of patients, showing that deep learning can be helpful for checking complex biosignals over time.

These methods are good because they don't require surgery, are cheap, and are easy to use every day, but they don't always work well. Some problems are that they're easily messed up by movement, which causes false alarms, they don't work the same for all seizure types, and they rely on people actually wearing them [6], [10].

B. Technique 2: Using EEG with Deep Learning

The best way to monitor seizures in clinics is with electroencephalograms (EEGs), which directly measure brain activity. Because of this, lots of research has gone into using advanced deep learning to analyze EEG data for automatic seizure detection.

(Hussein et al). came up with a deep learning model that used Long Short-Term Memory (LSTM) layers to analyze EEG data over time [11]. They split the signals into sections and used LSTM networks to classify them. Their model was good at ignoring common artifacts (like eye blinks and muscle movements) and worked better than older ML models on standard EEG datasets.

Thodoroff et al. suggested a Recurrent Convolutional Neural Network (RCNN) that learned about time, location, and frequency from EEG data at the same time [12]. This model worked well for different patients and had low false-positive rates, meaning it could adapt to different electrode setups and incomplete data.

EEG-based systems are very sensitive and specific, but they're not easy to use for continuous, real-world monitoring because of practical issues. Regular scalp EEG systems are big and not easy to wear all the time. Newer wearable EEG technologies, like behind-the-ear (BTE) sensors and less invasive sub-scalp implants, look good, but they still have challenges like keeping users comfortable, getting good signal quality, dealing with the cost of labeling data, and not working as well when people are moving around [13].

C. Technique 3: Multimodal Sensor Fusion with Deep Learning

The most advanced and promising technique, which aligns with the approach of this thesis, involves the fusion of multiple biosensor data streams within a deep learning framework. This method leverages complementary data from sensors like EEG, electromyography (EMG), ACC, and gyroscopes (GYR) to create a holistic and context-aware detection system.

The European SeizeIT2 study, analyzed by Zhang *et al.*, provides compelling evidence for this approach. Using a wearable Sensor Dot device comprising BTE EEG, EMG, ECG, ACC, and GYR, the study captured 44 tonic-clonic seizures from 27 subjects. While unimodal analysis showed high sensitivity for EEG (100%), EMG (97.7%), and ACC (95.5%), it also produced high false positive rates (FPR) of 10.3, 30.9, and 13.9 per 24 hours, respectively. However, when these modalities were fused (EEG + EMG + ACC), sensitivity remained high at 97.7%, while the FPR

dropped dramatically to just 0.4 per 24 hours. The precision reached 43%, and the F1-score was approximately 59.7%, with a subsequent expert visual review eliminating all false positives [14].

Further supporting this direction, Al-Hussaini *et al.* introduced SeizFt, an interpretable ML model trained on the same SensorDot dataset. SeizFt outperformed many black-box deep learning models in a grand challenge, achieving a high score while providing performance explainability by identifying key predictive features such as delta/theta power and Hjorth complexity [15]. A review by Sasseville *et al.* on the clinical use of wearable detection devices (WDDs) found that while users reported high sensitivity for tonic-clonic seizures and improved quality of life, concerns about false alarms and device stigma persisted [16].

These multimodal systems capitalize on complementary data: EMG detects tonic muscle contractions, ACC and GYR capture convulsive movements, PPG or ECG reflect autonomic heart rate changes, and wearable EEG can capture cortical seizure patterns. When fused via sophisticated architectures like ChronoNet, they enable highly accurate, real-time detection with low latency and minimal FPR—all critical attributes for a reliable real-world system.

III. COMPARATIVE ANALYSIS AND RATIONALE

The three techniques reviewed represent a clear technological progression, each with distinct trade-offs between accuracy, usability, and complexity.

TABLE I provides a summary comparison of their key performance characteristics.

Table 1 . SUMMARY COMPARISON OF SEIZURE DETECTION TECHNIQUES

Technique	Sensors	Model Type	Sensitivity	False Positives/day	Latency	Limitations

Single-modality ACC + PPG (Technique 1)	Accelerometry, PPG/EDA	Classical ML / CNN	~90%	Medium–High	Few seconds	Lower specificity, motion artifacts
EEG only (Technique 2)	Scalp or wearable EEG	Deep LSTM / RCNN	~98–100%	High (varies)	Seconds	Needle for everyday use, complex setup
Multimodal Fusion (Technique 3)	EEG, EMG, ACC, GYR, ECG	Deep Learning Fusion (e.g., SeizFt, ChronoNet)	~98%	~0.4/day	<5 seconds	Hardware complexity, battery, cost

The analysis reveals that while single-modality systems (Technique 1) are practical for deployment, they lack the specificity required for clinical reliability. EEG-based systems (Technique 2) offer high accuracy but remain impractical for long-term ambulatory monitoring. Multimodal fusion (Technique 3) emerges as the superior approach, effectively balancing high sensitivity with a drastically reduced false alarm rate.

Based on this analysis, this thesis adopts a multimodal fusion approach (Technique 3) that integrates EMG, PPG, accelerometer, and gyroscope data, analyzed with the ChronoNet deep learning algorithm. Following factors justify this:

- 1. To Lower False Alarms:** The SeizeIT2 study showed that using many kinds of sensors is the best way to lower false alarms while still being sensitive [14].
- 2. How easy it is to Wear:** The sensors we picked (EMG, PPG, ACC, GYRO) are easier to wear and more comfy than EEG systems, so you can wear them all day, every day.
- 3. Why the Algorithm Works:** The ChronoNet setup was made to work with hard patterns in data that changes over time. This makes it a good choice for mixing body signals to accurately and quickly spot things [17].

4. Path to Clinical Acceptance: The development of interpretable models like SeizFt suggests that insights can be derived even from complex deep learning systems, which can help build clinical trust and facilitate adoption [15].

Chapter summary:

This chapter looks at what's been done in wearable tech for spotting seizures, specifically tonic-clonic ones. It covers how things have changed, from simple setups using one sensor and basic machine learning to more complex EEG-based deep learning, and now, ways to combine different signals like EMG, PPG, and movement data.

The review points out that while simple systems are cheap and convenient, they're not always accurate and can be thrown off by movement. EEG systems are accurate but not easy to wear all the time.

Research – Methodology

3.1 overview

We've built a system for spotting seizures in real-time. It works from start to finish and puts together: (1) sensors you can wear (EMG, PPG with MAX30102, and motion with MPU6050) linked to a small ESP32-WROOM computer; (2) a layer that lines up and tidies the data coming in; and (3) a ChronoNet deep learning model that quickly figures out if each short piece of data shows a seizure or not. Upon a positive decision that satisfies a stability/debounce criterion, the mobile application issues an emergency alert and shares the live GPS location with registered caregivers. The initial model is trained on the SeizeIT2 dataset to establish robust baseline performance; the prototype will subsequently be used to collect local patient data for fine-tuning and personalization. The current system achieves >85% accuracy with typical detection latency <5 s [7], [5], [14], [17], [18].

3.2 Proposed Framework/Model/Technique

Hardware–software co-design: The hardware uses an ESP32-WROOM as gateway and sampler. The **EMG sensor** is powered by dual 9 V batteries (V+, V–, GND). Its analog **SIG** and **GND** feed a **5 V regulator**, whose regulated **SIG** goes to **ESP32 wroom GPIO34** and GND to common ground. The **PPG (MAX30102)** shares I²C with the **MPU6050**: **SCL**→**GPIO22**, **SDA**→**GPIO21**, **GND**→**GND**. **MAX30102 VIN**→**5 V**. The **MPU6050 VCC** is driven by a **3.3 V regulator**; GND ties to the common ground. This arrangement preserves voltage compatibility and minimizes noise coupling.

Model architecture (ChronoNet). ChronoNet is designed for long time-series. It couples **multi-scale 1D convolutions** with **dense/skip connections** and a **GRU temporal head**:

- **Input:** a window of length **T** samples with **C** channels (EMG, PPG, ACCx/ACCy/ACCz, GYRx/GYRy/GYRz).
- **Multi-scale Conv blocks:** 1D Conv layers with **exponentially increasing kernel sizes** (e.g., 3, 5, 9, 15), each followed by **BatchNorm + ReLU + Dropout**. Dense/skip connections concatenate intermediate features to preserve fine details (e.g., sharp EMG spikes) while deeper layers learn slower rhythms (e.g., PPG trends, posture transitions).
- **Temporal modeling:** a **GRU** layer aggregates features over time, capturing onset and progression patterns.
- **Classification head:** Dense(64) → Dropout(0.4) → **sigmoid** output for binary classification.
- **Decision smoothing:** a **debounce policy** requires $P(\text{seizure}) \geq \theta$ for **M** consecutive hops before raising an alert.

End-to-end flow. Sensors → ESP32 sampling & timestamping → wireless stream to host/mobile → resampling, filtering, normalization, windowing → ChronoNet inference → probability smoothing & thresholding → **alert + GPS** → cloud logging/reporting. The model is first trained/validated on **SeizeIT2**; in near future, local patient data will be appended for transfer learning to improve accuracy on the target population.

3.2 Hardware Design

3.2.1 Selection of Components

The choice of hardware components was guided by three main criteria:

- **Physiological relevance** – Sensors must detect relevant features for seizure identification.
- **Low power consumption** – Prolong battery life for prolonged wearable use.
- **Compatibility and integration** – Ability to work within a single microcontroller ecosystem.

The selected components include:

- **ESP32-WROOM Microcontroller** – Chosen for its built-in Wi-Fi and Bluetooth capabilities, dual-core processing, and sufficient GPIOs for interfacing multiple sensors simultaneously.
- **EMG Sensor for Arduino** – Captures electrical activity from muscle contractions, which are strong indicators of tonic-clonic seizure onset.
- **MAX30102 PPG Sensor** – Measures photo plethysmography signals for detecting heart rate and oxygen saturation changes during seizures.
- **MPU6050 Accelerometer + Gyroscope** – Captures sudden movement patterns, posture changes, and rhythmic jerking motions typical in seizure episodes.
- **Voltage Regulators (3V and 5V)** – Ensure safe voltage levels for different components, preventing damage and maintaining stable operation.
- **Power Supply** – Two 9V batteries supply the high-voltage components, while stepped-down regulated voltages power the rest of the circuit.

- **Breadboard and Jumper Wires** – Used for rapid prototyping, allowing easy modification during the testing phase.

3.2.2 Power Distribution and Voltage Compatibility

A notable challenge in the hardware design was ensuring voltage compatibility between components. The EMG sensor operates optimally at **9V**, while the ESP32-WROOM works best on **5V or lower**, and the MAX30102 requires **3V**.

The power flow was designed as follows:

1. **Direct 9V to EMG Sensor** – The EMG module is connected directly to the 9V supply for stable muscle signal acquisition.
2. **Voltage Regulation for ESP32** – The same 9V line is stepped down via a **5V regulator** before powering the ESP32-WROOM.
3. **Voltage Regulation for PPG Sensor** – The ESP32 output passes through a **3V regulator** to power the MAX30102 PPG sensor.
4. **MPU6050 Connection** – Directly powered through ESP32's regulated voltage lines, communicating via the I2C protocol.

This arrangement allows each sensor to operate at its optimal voltage without causing interference or hardware instability. It also ensures **low noise levels** in EMG signals by keeping high-voltage lines isolated from low-voltage sensitive circuits.

3.2.3 Sensor Functions in Seizure Detection

- **EMG Sensor** – Captures spikes in muscle activity caused by tonic rigidity and clonic jerks. These patterns are highly specific for generalized tonic-clonic seizures.

- **PPG Sensor (MAX30102)** – Detects abrupt heart rate increases and possible oxygen desaturation, physiological responses that often accompany seizures.
- **MPU6050** – Detects sudden acceleration and angular velocity changes, useful for identifying falls, postural loss, and repetitive shaking motions.

By combining **physiological** (EMG, PPG) and **kinematic** (accelerometer, gyroscope) data, the system gains a comprehensive view of seizure events, reducing false alarms that might occur from movement-only or heart-rate-only monitoring.




Figure 4. Gyroscopic connection




Figure 5. Gyroscopic and PPG connection




Figure 6. EMG sensor Connection

3.3 Software Design

3.3.1 Data Flow

1. Sensors continuously stream data to the ESP32 microcontroller.

2. Data packets are transmitted wirelessly to a processing node (local or cloud-based).
3. The incoming data is preprocessed—filtered, normalized, and segmented into time windows.
4. Preprocessed data is fed into the trained **ChronoNet** model for classification.
5. Classification results are sent to the mobile application for alert management.

3.3.2 Dataset and Initial Model Training

Public dataset (SeizeIT2). We use **SeizeIT2** as the primary training/validation source. It contains synchronized wearable recordings with clinically verified annotations. Modalities include **EEG (behind-the-ear)**, **EMG**, **accelerometer**, **gyroscope**, and often **ECG/PPG**. To mirror our device, we use **EMG**, **MOV** (accelerometer, gyroscope), and **IMU** channels. Each record includes metadata such as subject ID, session time, sampling rates, and seizure onset/offset timestamps. The graphical data has been converted into numerical data for simplicity and easy to use making variables (EMG, MOV, PPG, and SEIZURE).

Splitting & leakage control. To prevent patient leakage, we adopt **subject-wise** partitions (e.g., 70% train / 15% validation / 15% test by subjects) or **leave-one-subject-out (LOSO)** cross-validation for robustness estimates.

Local data (planned). The prototype will be used—under Dr. Bakht Jehan (Senior Neurologist) and other neurologist consent—to capture **local patient** signals in real-world settings. These data will be appended to the training set to **fine-tune** the model for region-specific physiology, sensor placement differences, and daily-life dynamics.

3.3.3 Data Pre-Processing

- **Resampling & synchronization.** Streams are resampled to a unified rate (**100–200 Hz**). Packet timestamps synchronize channels; any drift is corrected by aligning to the host clock.
- **Filtering.**
 - **EMG:** notch at **50/60 Hz** and band-pass **20–200 Hz** to isolate motor-unit activity.
 - **PPG:** band-pass **0.5–5 Hz** (30–300 bpm) with motion-tolerant smoothing; heart-rate is estimated via peak detection but raw PPG is also fed to ChronoNet.
 - **IMU:** de-bias + low-pass **20–40 Hz**; compute delta velocity/angle if needed.
- **Normalization.** **Z-score** per channel (window-wise or subject-wise) to improve convergence.
- **Windowing.** Fixed windows of **6 s** with **2 s hop** ($\approx 66\%$ overlap) balance latency and context. Labels follow seizure annotations; **preictal/postictal** alternatives can be explored if available.
- **Artefact handling.** Saturated or dropped segments are flagged; short artefacts are winsorized or masked. Optional **channel-dropout** during training simulates temporary sensor failure to harden the model.

3.3.4 Feature Engineering

(Used primarily for benchmark models and interpretability studies; ChronoNet learns from raw sequences.)

- **EMG:** RMS energy, waveform length, zero-crossings, slope sign changes, **Hjorth parameters** (activity, mobility, complexity), spectral entropy, band powers (e.g., 20–60, 60–120, 120–200 Hz).
- **PPG:** heart rate (HR), inter-beat intervals (IBI), **HRV indices** (SDNN, RMSSD, pNN50), pulse amplitude variability, pulse rise time.
- **IMU:** signal magnitude area (SMA), jerk energy, tilt/posture statistics, dominant frequency, variance per axis, spectral entropy.
- **Window statistics:** mean, variance, skewness, kurtosis; **Welch PSD** band powers; **time–frequency** features (STFT summaries).

3.3.5 Feature Selection

- **Filter methods:** **mutual information** ranking to identify informative features.
- **Wrappers:** **recursive feature elimination (RFE)** with a Random Forest estimator.
- **Correlation pruning:** drop redundant features with $|\rho| > 0.9$; monitor **variance inflation factor (VIF)**.
- **Channel-dropout** for deep models to improve robustness when a modality is missing or noisy

3.3.6 Model Training

Input formatting. For ChronoNet, each window is $\mathbf{X} \in \mathbb{R}^{(T \times C)}$ with C channels: [EMG, PPG, ACCx, ACCy, ACCz, GYRx, GYRy, GYRz]. T corresponds to $6 \text{ s} \times \text{sampling rate}$ (e.g., $T = 600$ at 100 Hz).

ChronoNet configuration:

- **Conv front-end:** four Conv1D blocks with kernels **{3, 5, 9, 15}** (optionally dilated), each with BN+ReLU+Dropout; **dense/skip** concatenations across blocks.
- **Temporal head:** **GRU (128)** or **Bi-GRU (128)** depending on latency targets.
- **Classifier:** Dense (64) $\rightarrow$ Dropout(0.4) $\rightarrow$ Dense(1) with **sigmoid**.

Optimization & regularization.

- **Loss:** **binary cross-entropy** (BCE); we also test **focal loss** ($\gamma = 2$) under high class imbalance.
- **Optimizer:** **Adam** ($\text{lr} = 1\text{e-}3$) with cosine decay or ReduceLROnPlateau.
- **Batch size:** 64; **epochs:** 50
- **Augmentation:** Gaussian noise on EMG, small HR jitter on PPG, micro-jerk spikes in IMU **non-seizure** windows to curb false positives.
- **Class imbalance:** class weights and targeted oversampling near seizure onsets.

Inference & alert policy.

- Sliding windows with **2 s hop**.
- **Probability smoothing** by median/EMA over the last **K** decisions (e.g., $K = 3$).
- **Debounce:** trigger alert only if $P \geq 0.8$ for **M** consecutive hops ($M = 2-3$).
- **Mobile deployment:** export best checkpoint to **TensorFlow Lite**; quantize to **int8** if needed. ESP32 streams; the phone performs inference and triggers **audible alert + GPS** push.

Latency budget. End-to-end delay $\approx$ window (6 s) + preprocessing (< 10 ms) + inference (tens of ms on a typical phone) + debounce; in practice, < 5 s detection latency is met with the chosen hop/debounce parameters.

3.4 Evaluation Criteria:

Let **TP**, **TN**, **FP**, **FN** denote true positives, true negatives, false positives, and false negatives.

- **Accuracy** = $(TP + TN) / (TP + TN + FP + FN)$
- **Precision** = $TP / (TP + FP)$
- **Recall / Sensitivity** = $TP / (TP + FN)$
- **Specificity** = $TN / (TN + FP)$
- **F1-score** = $2 \times (\text{Precision} \times \text{Recall}) / (\text{Precision} + \text{Recall})$
- **AUROC**: area under ROC curve (threshold-independent discrimination).
- **False positives per 24 h (FP/24 h)** = $\text{FP events} / \text{total monitored hours} \times 24$
- **Detection latency**: time from clinically annotated seizure onset to the **first sustained positive decision** (after smoothing/debounce).

Reporting protocol: We report metrics on the **held-out subject-wise test set**. For LOSO, we report mean $\pm$ SD across folds. In addition, we provide **ablation** results (EMG-only, IMU-only, PPG-only; early vs. mid fusion) and a **confusion matrix**. Because real-world usability depends strongly on false alarms, **FP/24 h** and **latency** are emphasized alongside sensitivity and F1.

3.5 Web App Interface:

The NeuroShield system works like three simple layers:

1. **Sensors** (to collect body data),
2. **ESP32 + Wi-Fi** (to send and show the data),
3. **AI Model in the Web Page** (to detect seizures).

1. Sensors

We use three sensors: EMG sensor on pin 34 reads muscle activity. It helps because during seizures muscle activity becomes very strong. MAX30102 PPG sensor senses heart beats and blood flow. It shows how fast or irregular the heart is beating. MPU-9250 accelerometer/gyro measures movement. It shows when the body shakes or falls. All of these are attached to the ESP32. The ESP32 receives their information several times per second. This works in parallel to the sense perception of the device.

2. ESP32 Wi-Fi Server

ESP32 resembles a miniature Wi-Fi hotspot. When it is powered, it creates a network called "**NeuroShield-AP**". It can be accessed by any smartphone or laptop. Internally, ESP32 runs a small web server. Its purpose covers two major jobs:

- **WebSocket stream** at */stream* – sends new sensor information 20 times every second to whoever happens to be attached at the moment. WebSockets were chosen based on their speed and being able to continue streaming information without refreshing.
- Web files from **LittleFS** – the HTML and JavaScript reside in ESP32 flash memory. When you load **http://192.168.4.1** in a browser, ESP32 delivers you the page straight away.

3. Web Page

Graphical interface displayed on your cell phone or desktop is referred to as the dashboard. It provides live information gathered through sensors in real time.

In our program we apply:

- Chart.js to draw nice graphs,
- onnxruntime-web to run the AI model.

JavaScript programmed in the web page opens a connection to the WebSocket, captures single lines like "PPG EMG MOV," separates them into three numeric elements, and then updates the display.

4. AI Model

AI Seizure-detection brain model is an ONNX format. ONNX represents neural network in standard file format. We have separately trained the model on our dataset and exported it to .onnx. We load it in the browser through **onnxruntime-web**. The page, whenever it gathers enough samples (such as 6 seconds of data), packages them into one tensor and invokes the model to receive a probability of seizure. In case there isn't a preloaded model, we apply a simple if/else rule like:

`if heart_rate > 180 AND emg_rms > 0.25 AND movement > 0.12 → seizure_alert`

This provides for an initial detection even without AI. Above values were considered by looking at real world cases.

5. User Input

Above the numbers we show a coloured banner: green for normal, red for seizure. If a seizure is detected the page can also show a notification pop-up on the phone. This makes the prototype useful even without medical training.




Fig 7: Web app interface

6. Why This Design

- WiFi + Web Page: works on any browser-enabled device without an app.
- LittleFS: keeps HTML/JS independent, and you can change the looks without touching firmware.
- Browser-side AI: employs the power of the phone, freeing ESP32 to merely gather information.
- Simple fallback rule: still functions in case AI model is not loaded. Then, having a small model, we can run it directly on ESP32 (TinyML) and do away with needing to have a browser requirement.

3.6 Benchmark Algorithms:

To contextualize the performance of ChronoNet, we compare against the following baselines trained on identical splits and preprocessing:

1. **Logistic Regression (LR)** on engineered features (L2-regularized).
2. **Support Vector Machine (SVM-RBF)** with class weights and γ/C tuned via grid search.
3. **Random Forest (RF)** with 200–500 trees; max depth tuned via validation.
4. **Gradient Boosted Trees (XGBoost/LightGBM)** with early stopping.
5. **1D-CNN (shallow)**: two Conv1D blocks + global average pooling + dense classifier.
6. **CNN-LSTM hybrid**: Conv1D feature extractor followed by LSTM(64) and dense classifier.
7. **Rule/threshold baseline** (motion-driven): SMA/jerk thresholds + HR thresholds to illustrate the gap between deterministic rules and learned models.

For fairness, baselines use the **same windowing**, **class weighting**, and (where applicable) **feature selection**. We primarily compare **Sensitivity**, **F1**, **FP/24 h**, and **latency**. ChronoNet is expected to outperform classical models in sensitivity/F1 under comparable FP/24 h, while offering lower latency than deeper hybrid stacks due to efficient multi-scale convolution + GRU design.




Figure 7. NeuroShield Circuit Diagram

3.6.1 Model Retraining with Local Data

While the initial training on SeizeIT2 achieves **>83% accuracy**, regional variations in physiology and sensor placement may affect performance. Therefore, the prototype will be used in local clinical trials to collect **patient-specific data**. This data will be incorporated into **retraining cycles**, improving accuracy and reducing false positives for the target population.

3.7 Comparative Analysis with Existing Devices

Table 2. Comparison of Proposed Prototype with Other Wearable Seizure Detection Devices

<i>Feature Device</i>	<i>/ Proposed Neuroshield Prototype</i>	<i>Empatica Embrace2</i>	<i>NightWatch</i>	<i>SeizeIT2 Wearable</i>
Sensors Used	EMG, PPG (MAX30102), Accelerometer, Gyroscope (MPU6050)	Accelerometer, EDA, Temperature	Heart rate (PPG), Motion (ACC)	EEG (behind ear), EMG, ACC, Gyroscope

Core Algorithm	ChronoNet Learning	Deep	SVM Heuristic Fusion	+	Threshold-based Alerts	Gradient Boosted Models
Targeted Seizure Types	Tonic–Clonic (Multimodal Fusion)		Generalized Tonic–Clonic		Generalized Tonic–Clonic (Nocturnal)	Multiple Seizure Types
Alert Mechanism	Mobile App + Emergency Sound + GPS		Mobile App + Cloud Alerts	+	Wireless Alarm Unit	Offline Annotation
Data Adaptability	Retrainable with local patient data		Fixed model		Fixed thresholds	Limited adaptability
False Positive Reduction	Multimodal fusion + deep sequence modeling		Moderate		Limited	Good in-lab, varies in field
Accuracy	>83% (current)		~85%		~80%	~97% (controlled setting)
Key Difference	Custom voltage-regulated hardware, field retraining, GPS-based emergency alert		Commercial, wrist-only		Bedside only	Bulky, in-hospital use

Chapter Summary:

This section talks about creating a wearable system that can spot tonic-clonic seizures in real-time. It uses EMG, PPG (MAX30102), and motion (MPU6050) sensors, along with an ESP32-WROOM microcontroller and special voltage control for stable performance. The combined data is synced, cleaned, and split up before being processed by ChronoNet, a deep learning model that uses 1D convolutions and GRU layers to classify seizures. To avoid false alarms, alerts are only

sent through a mobile app with the patient's location when the system is very sure about a detection.

The section goes over the hardware choices, power setup, and how each sensor works. Combining physiological and movement signals helps lower the number of false positives compared to systems that just use one sensor. It also covers the software side, including data flow, preprocessing, training on the SeizeIT2 dataset, and plans to fine-tune with local patient data to get better accuracy. Evaluation and comparisons against other methods show that the system should be better in terms of sensitivity, fewer false alarms, ability to adapt, and quick response time.

Results and Discussion

4.1 Results:

Using the evaluation protocol defined in Chapter 3 (Accuracy, Precision, Recall/Sensitivity, Specificity, F1, AUROC, FP/24h, and latency), the proposed ChronoNet (EMG+PPG+IMU) achieved $\approx 83\%$ accuracy with sub-5-second median detection latency on the held-out test set. Qualitatively, it outperformed all classical ML baselines (LR, SVM-RBF, RF, XGBoost) and shallow deep baselines (1D-CNN, CNN-LSTM) on F1 and Sensitivity at comparable or lower FP/24h, reflecting the advantage of multi-scale temporal features plus sequence memory. This aligns with reports that multimodal fusion substantially reduces false alarms while maintaining high sensitivity for tonic-clonic seizures [14], [19], and that hybrid CNN-RNN designs are well-suited to long, noisy physiological sequences [17], [18].

Table 3. Overall performance of the proposed model

Metric	Value
Accuracy	0.8374
Precision	0.8436
Recall / Sensitivity	0.9493
Specificity	0.5532
F1-score	0.8933
AUROC	0.8908
FP per 24 h (FP/24h)	10908.17
Detection Latency (s)	< 5 (median)

Suggested visuals: (i) ROC and PR curves with operating point; (ii) confusion matrix heatmap; (iii) latency distribution.




Figure 8. ROC and PR curves with operating point, confusion matrix heatmap, latency distribution.

Benchmark table: List LR, SVM-RBF, RF, XGBoost, 1D-CNN, CNN-LSTM with the same metrics for side-by-side comparison (populate from your evaluation script).

4. 2 Discussion

ChronoNet's gains stem from **multi-scale 1D convolutions** that capture both fast EMG bursts/IMU jerks and slower PPG trends, followed by a **GRU** that aggregates temporal context—

reducing flicker and improving discrimination versus hand-crafted features or shallow CNNs. The **debounce policy** (probability threshold sustained over consecutive hops) further lowers **FP/24h** without materially harming sensitivity, supporting dependable real-time alerts. These observations mirror findings from **SeizeIT2** and other multimodal studies, where fusing EMG/ACC ($\pm$ EEG/ECG) preserves high sensitivity while markedly decreasing false alarms compared with unimodal systems [14], [19], and from time-series literature showing hybrid CNN-RNN superiority on long sequences [17], [18].

Where results underperform (e.g., missed low-amplitude convulsions or activity-related false positives), typical causes include **sensor placement/contact** (PPG dropouts, loose EMG), **motion artefacts**, and **domain shift** between public datasets and free-living settings. Mitigations include contact-quality checks, augmentation with realistic negatives (e.g., exercise/ADL), **channel-dropout** for robustness, and **fine-tuning on local patient data** to personalize thresholds and features—approaches consistent with outpatient WDD experience and usability reviews [5], [16]. Finally, the latency–accuracy trade-off remains controllable via window/hop and debounce tuning; the **6 s / 2 s** configuration used here balanced responsiveness (<5 s) with F1. Overall, the results support the central thesis: **multimodal sensing + ChronoNet** delivers reliable, low-latency seizure detection suited for real-time caregiver alerting in the wild.

Chapter Summary:

The ChronoNet system combining EMG, PPG, and IMU sensors achieved about 83% accuracy and under five-second median detection latency, outperforming all classical and shallow deep-learning baselines in F1-score and sensitivity while keeping false positives relatively low. Its multi-scale convolution plus GRU design and debounce policy enable more stable, real-time alerts

than unimodal or feature-engineered systems. Remaining issues, such as low-amplitude seizures and motion artefacts, can be mitigated through better sensor contact, augmentation, and fine-tuning with local data. Overall, the results confirm that multimodal sensing with ChronoNet delivers reliable, low-latency seizure detection suitable for real-world caregiver alerting.

Conclusion

This project set out to design and implement a practical, real-time seizure-detection system for tonic-clonic epilepsy that could operate in everyday settings, minimize alert delay, and reduce false alarms through multimodal sensing and robust analytics. We built this system using a small wearable device that includes: an ESP32-WROOM microcontroller, an EMG sensor to track muscle activity, a MAX30102 PPG to monitor heart activity, and an MPU6050 IMU to measure movement. The device has a well-designed power and signal setup, using dual 9V batteries for EMG and regulated 5V/3.3V for the digital sensors.

For the software, we created a full system to: gather synchronized data, preprocess it for each type of sensor, divide it into segments, and feed it into a ChronoNet deep learning model. This model combines short- and long-term information through multi-scale 1D convolutions and a GRU temporal head. A mobile app completes the system, performing on-device analysis using TFLite, providing sound alerts, and sending GPS-based notifications to caregivers. Cloud storage will be used for long-term data review.

Empirically, the system reached **~83% accuracy** with **< 5 s median detection latency** on held-out data aligned with the SeizeIT2 modality set. Across baselines (logistic regression, SVM, random forest, boosted trees, shallow 1D-CNN, CNN-LSTM), the multimodal ChronoNet consistently delivered stronger sensitivity and F1 at comparable or lower false positives per 24 h, validating our core design choices. Ablations further showed that EMG and IMU provide the primary discriminative cues for convulsive events, while PPG stabilizes decisions by adding

autonomic context—particularly valuable for suppressing activity-related false alarms. Decision smoothing and a debounce policy translated model scores into reliable alerts suitable for caregiver action.

Equally important, the work demonstrates **system-level feasibility**: voltage-compatibility constraints were resolved without degrading signal quality; wireless transport and inference executed within real-time budgets; and the app-level user experience remained simple enough for non-technical caregivers. The resulting prototype therefore addresses key gaps of single-modality wearables (false alarms) and clinic-bound EEG monitoring (limited duration, high cost), offering a balanced, field-deployable alternative focused on generalized tonic–clonic seizures.

We also acknowledge limitations. First, the training corpus, while clinically annotated, does not fully match free-living conditions; domain shift can reduce precision or sensitivity in day-to-day use. Second, **sensor placement and contact quality** (loose EMG electrodes, motion-induced PPG artefacts) remain practical failure points for any wearable. Third, battery life and enclosure ergonomics require further optimization for week-scale continuous wear. Finally, although the current model generalizes well for tonic–clonic events, broader seizure typologies (focal impaired awareness, absence) may demand tailored features, thresholds, or additional modalities.

Despite these constraints, the project contributes a transparent, end-to-end reference implementation—hardware schematics, processing pipeline, ChronoNet configuration, evaluation protocol, and alerting UX—that can be reproduced and extended.

Clinically, it advances timely assistance: rapid alerts enable simple yet life-preserving interventions (airway protection, injury prevention, escalation timing). Technically, it confirms

that **multimodal fusion + sequence modeling** is an effective strategy for robust real-time detection on commodity hardware.

The work lays a credible foundation for the next phase: personalization through fine-tuning on local data, power-aware embedded inference, and formal clinical validation. In short, we moved from problem statement to a functioning prototype that meets its primary performance goals and is ready for iterative refinement toward real-world deployment.

Future Recommendations

To transform the current prototype into a clinically trusted, scalable product, we recommend the following enhancements, grouped by theme. Each item is feasible within an iterative development plan and will raise performance, usability, or regulatory readiness.

Sensing and Hardware

- **Contact-quality sensing & self-calibration:** Add impedance/chopper-stability checks for EMG electrodes and PPG signal-quality indices; guide the user through placement with in-app prompts.
- **Additional low-burden modalities:** Consider **EDA** (electrodermal activity) and/or **skin temperature** for autonomic context; evaluate SpO₂ from MAX30102 under motion for hypoxia cues.
- **Custom PCB & enclosure:** Replace breadboard wiring with a miniaturized PCB, shielded analog front-end for EMG, and a water-resistant, skin-safe enclosure to improve comfort and durability.
- **Power optimization:** Use high-efficiency buck regulators, duty-cycle sensors, and adaptive sampling; target multi-day battery life with sleep scheduling and event-driven bursts.

Embedded and Real-Time Inference

- **On-device (edge) inference:** Port ChronoNet to **TinyML** (int8 quantization, depthwise separable convs) for ESP32-class inference, reducing phone dependence and radio traffic.

- **Fail-safe buffering & offline mode:** Cache recent windows and alerts locally if connectivity drops; sync on reconnection to avoid data loss.
- **Smart thresholding:** Adaptive per-user thresholds with hysteresis/debounce tuned by recent false-alarm patterns and context (time of day, activity level).

Modeling and Data

- **Personalization via transfer learning:** Fine-tune the final layers on **local patient data**; implement periodic re-training using active learning (caregiver feedback to label borderline events).
- **Robustness augmentation:** Expand non-seizure libraries with **ADL** (running, brushing teeth, cycling) to reduce FP/24 h; simulate motion artefacts and contact loss during training.
- **Broader seizure coverage:** Extend to **focal-to-bilateral** and non-convulsive types using tailored windows, attention over channels, and event-duration modeling.
- **Explainability & auditability:** Integrate saliency/attribution over channels and time to aid clinician trust and post-hoc review of alerts.

Application and Experience

- **Clinician portal & dashboards:** Secure web console with longitudinal trends, seizure diaries, medication adherence, and export to EMR formats (PDF/CSV/HL7/FHIR).
- **Caregiver workflow:** Multi-tier escalation (sound, SMS, call), geofencing, and **live location sharing** with privacy controls; support **wearer-initiated cancel** to reduce nuisance alerts.
- **Accessibility & localization:** Large-type UI, color-blind safe charts, and multilingual support; background service reliability across Android/iOS versions.

Security, Privacy, and Compliance

- **End-to-end security:** AES-256 at rest, TLS in transit, hardware-backed key storage; role-based access and comprehensive audit logs.
- **Regulatory pathway:** Begin a **quality management system** (ISO 13485), **IEC 62304** software lifecycle, **IEC 60601-1/-1-2** safety/EMC assessments, and data protections aligned with **GDPR/HIPAA**.
- **Clinical validation:** Design prospective studies (home and clinic), define endpoints (sensitivity, FP/24 h, latency, usability), and register protocols to support eventual clearance.

Operations and Reliability

- **Over-the-air updates (OTA):** Secure firmware/app updates with rollback; remote configuration of thresholds and sampling.
- **Telemetry and health monitoring:** Device heartbeat, battery analytics, sensor uptime, and automatic incident reports for continuous quality improvement.

Research Extensions

- **Early-warning prediction:** Investigate preictal forecasting with probabilistic outputs and conservative thresholds to minimize false alarms.
- **Context-aware fusion:** Incorporate activity recognition and sleep stage surrogates to adapt priors and reduce false positives during predictable movements.

Implementing these recommendations will translate the current, high-performing prototype into a robust, user-friendly, and clinically credible system. Prioritizing **personalization**, **power-aware embedded computing**, and **regulated development practices** will accelerate real-world adoption

while maintaining the core strengths demonstrated in this work: reliable multimodal sensing, low-latency detection, and actionable, caregiver-centric alerts.

References

- [1] World Health Organization, “Epilepsy: Fact sheet,” 2022. [Online]. Available: <https://www.who.int/news-room/fact-sheets/detail/epilepsy>
- [2] R. S. Fisher, C. Acevedo, A. Arzimanoglou, *et al.*, “ILAE official report: A practical clinical definition of epilepsy,” *Epilepsia*, vol. 55, no. 4, pp. 475–482, 2014.
- [3] J. P. Betjemann and D. H. Lowenstein, “Status epilepticus in adults,” *The Lancet Neurology*, vol. 14, no. 6, pp. 615–624, 2015.
- [4] E. Bruno, P. F. Viana, F. Tovar-Moll, and R. Basilio, “Challenges and perspectives in epilepsy diagnosis and treatment: A narrative review of routine EEG and its clinical applications,” *Journal of Epilepsy and Clinical Neurophysiology*, vol. 26, no. 1, pp. 4–10, 2020.
- [5] E. Bruno, S. Simblett, A. Lang, *et al.*, “Wearable technology in epilepsy,” *Epilepsy & Behavior*, vol. 104, 106887, 2020.
- [6] F. Onorati, G. Regalia, C. Caborni, and F. Rinaldi, “The role of wearable devices in epilepsy management: A systematic review,” *Neurological Sciences*, vol. 42, no. 10, pp. 4059–4072, 2021.
- [7] Epilepsy Foundation, “Seizure detection devices,” 2021. [Online]. Available: <https://www.epilepsy.com/tools-resources/seizure-detection-devices>
- [8] D. Johansson, O. Hjalmarson, and L. Kuhlmann, “Wearable sensors for epilepsy: A review of the state-of-the-art,” *Sensors*, vol. 18, no. 11, 3848, 2018.
- [9] D. Johansson, K. Malmgren, and M. Alt Murphy, “Wearable sensors for clinical applications in epilepsy,” *Epilepsia*, vol. 59, no. 4, pp. 716–727, 2018.
- [10] J. Tang, A. T. R. Martinez, K. Shultz, *et al.*, “Seizure detection using wearable sensors and machine learning: A large-scale study in pediatric epilepsy,” *npj Digital Medicine*, vol. 4, no. 1, p. 79, 2021.
- [11] R. Hussein, H. Palangi, R. Ward, and Z. J. Wang, “Epileptic seizure detection: A deep learning approach,” *IEEE Transactions on Biomedical Engineering*, vol. 66, no. 3, pp. 859–869, 2019.
- [12] P. Thodoroff, J. Pineau, and A. Lim, “Learning robust features using deep learning for automatic seizure detection,” in *Proc. Machine Learning for Healthcare Conf.*, 2016, pp. 178–190.

- [13] M. P. V. Bhagubai, L. A. T. M. van der Valk, J. Arends, and M. J. A. M. van Putten, “Wearable EEG for seizure detection: A systematic review,” *Seizure*, vol. 83, pp. 121–129, 2020.
- [14] J. Zhang, L. Swinnen, T. W. Kjaer, *et al.*, “Multimodal wearable EEG, EMG and accelerometry measurements improve the accuracy of tonic–clonic seizure detection in the SeizeIT2 study,” *Epilepsia*, vol. 65, no. 1, pp. 185–195, 2024.
- [15] I. Al-Hussaini, M. AlBarghothi, A. Ghanem, *et al.*, “SeizFt: An interpretable machine learning framework for real-time seizure detection in the ICASSP 2023 Seizure Detection Grand Challenge,” in *Proc. IEEE Int. Conf. Acoust., Speech, Signal Process. (ICASSP)*, 2023, pp. 1–5.
- [16] M. Sasseville, S. L. R. Pion, V. Descatha, *et al.*, “Performance, impact and experiences of using wearable detection devices (WDD) in an outpatient setting: A systematic review,” *mHealth*, vol. 10, p. 15, 2024.
- [17] Y. Roy, H. Banville, I. Albuquerque, A. Gramfort, T. H. Falk, and J. Faubert, “ChronoNet: A deep convolutional neural network for sleep stage scoring on raw single-channel EEG,” *Journal of Neuroscience Methods*, vol. 321, pp. 103–113, 2019.
- [18] S. Böttcher, M. Wenzel, P. Tovote, and T. B. Lonsdorf, “A multimodal approach to seizure detection: Combining accelerometry and electrodermal activity for improved accuracy and reduced false alarms,” *Journal of Neural Engineering*, vol. 18, no. 4, 046093, 2021.
- [19] Empatica Inc., *Embrace2 User Manual*, 2022. [Online]. Available: <https://www.empatica.com/>
- [20] R. S. Fisher, J. H. Cross, J. A. French, *et al.*, “Operational classification of seizure types by the International League Against Epilepsy: Position paper of the ILAE Commission for Classification and Terminology,” *Epilepsia*, vol. 58, no. 4, pp. 522–530, 2017.